

SCAN•MATE™

BY ARMOR DENTAL

YOUR INTRAORAL SCANNING PARTNER

Comfortable retraction, control, and efficiency that produces quality images




Improving digital work flow for intraoral scanning

Error-free scans from the beginning

Latex-free and fully autoclavable instrument up to 60 cycles

Professional isolation, retraction instrument for intraoral scanning and any retraction needs

- Ultimate retraction instrument for all dental procedures or exams.
- Non-reflective device for enhanced scanning for single restoration to full arches.
- Ideal for all intraoral scanning needs – crowns, bridges, implants, and orthodontics.
- Comfortable retraction for scanners, lasers, intraoral cameras, and air polishers.
- Durable ergonomic handle to improve control and quality of image.
- Soft, bendable tip provides smooth retraction and easier access for hard to reach areas.



“Best scan I have ever taken with a retraction tool.”

DR. CARL BOSCKETTI



Scan•Mate is changing everything for digital work flow

- Reduces scanning time for dental team and results in cleaner scans.
- Improves patient experience and first impression chairside.
- More efficient, results in high quality scans.
- No more tissue issues when scanning.
- Helps the dental team improve their scanning skills, freeing up dentists' time to do other procedures.

SCAN•MATE™


BY ARMOR DENTAL

Comfortable for clinician and patient

- Scan•Mate can be used for scanning and all retraction needs.
- No more fingers or mirrors or other isolation products needed.
- Patented bendable neck and ergonomic handle make it easy to scan images with comfort – bite image, 2nd posterior molar, lower distal, unruly tongue, tight lips.
- Reduces clean-up time after scanning.
- Latex-free.

“Scan•Mate allows me to quickly and clearly visualize the area to be scanned, helping to insure accurate, crisp, and clean images.”

DR. SCOTT GANZ




To Autoclave Scan•Mate

- Steam autoclave at 250F to 273°F (121°C to 134°C) for at least 20 minutes as per ADA guidelines. (Time may vary depending on size of load and autoclave.)
- Autoclavable 60 cycles and then replace after 1 month.
- Unique bendable handle and soft tip improves the patient experience.
- Scan•Mate's adjustable neck will eventually weaken after a month.
- Do not keep over the recommended autoclave limit to ensure product effectiveness.

4 units per box/\$49.95, Item No. 50012
Order via your distributors.
Contact: info@armor-dental.com.

Manufactured by:
Armor Dental Corporation, Hingham, MA 02043
Patent No. 9,078,723. Made in USA.
www.armor-dental.com



Get your scan right
the first time.
You're welcome!

ARMORDENTAL